

Disposable Grasping Forceps (Retrieval net)

A Xinwell Member Company


AteTec

Intended use

It is used to retrieval foreign bodies and polyps from digestive tract under endoscope.

Products Features

- The handle and net bag can be rotated synchronously, the control is closer and more efficient.
- The mesh net design reduces the risk of specimen slippage and loss.
- The retrieval net adopts imported nylon conforming material, the contact surface is smoother, and it is easier to pass the physiological bending and narrowness of the digestive tract.
- The deep elastic bag design can effectively filter the liquid that the human body can grab naturally while improving the grasping and control of the specimen.



Model	Diameter of outer tube (mm)	Length (mm)	Opening width (mm)
ATE-YWQ-TY-18*1600*20	1.8	1600	20
ATE-YWQ-TY-23*1600*20	2.3	1600	20
ATE-YWQ-TY-23*1600*25	2.3	1600	25
ATE-YWQ-TY-23*1600*30	2.3	1600	30
ATE-YWQ-TY-23*1600*40	2.3	1600	40
ATE-YWQ-TY-23*2300*20	2.3	2300	20
ATE-YWQ-TY-23*2300*25	2.3	2300	25
ATE-YWQ-TY-23*2300*30	2.3	2300	30
ATE-YWQ-TY-23*2300*40	2.3	2300	40

JIANGSU ATE MEDICAL TECHNOLOGY CO.,LTD

ADD: No.8,Lanxiang Road,Wujin Economic Development Zone
213000 Jiangsu China

E-mail: atemedical@xinhaigroup.com Web: www.ate-medical.com

+86-519 8528 8988

+86-519 8528 8060

A Xinwell Member Company

AteTec