

With alligator jaws



without spike

Disposable Biopsy Forceps

Easy-Effective-Exactly

A Xinwell Member Company









AteTec

Intended use

The disposable biopsy forceps works with flexible endoscopes, passing through the endoscope channel into the human body cavity to take the living tissues for pathology analysis.

Products Features


- The jaws are made of imported special high-strength medical stainless steel, with sharp bite and less bleeding.
- The unique process is trimmed, the jaws are highly consistent, and the tissue is obtained more fully.
- Ergonomic handle design, comfortable and convenient operation, stronger controllability.
- The specifications and models are rich and varied, coated/uncoated; oval cup/alligator cup; with needle/without needle, meeting various clinical needs.

	Model	Diameter of jaws(mm)	Length(mm)
	Uncoated,without spike		
	ATE-QYQ-A-18×1050	1.8	1050
	ATE-QYQ-A-18×1200	1.8	1200
	ATE-QYQ-A-18×1600	1.8	1600
	ATE-QYQ-A-23×1600	2.3	1600
	ATE-QYQ-A-23×1800	2.3	1800
	ATE-QYQ-A-23×2300	2.3	2300
	Uncoated,with spike		
	ATE-QYQ-B-23×1600	2.3	1600
	ATE-QYQ-B-23×1800	2.3	1800
	ATE-QYQ-B-23×2300	2.3	2300
	Coated, without spike		
	ATE-QYQ-C-18×1050	1.8	1050
	ATE-QYQ-C-18×1200	1.8	1200
	ATE-QYQ-C-18×1600	1.8	1600
	ATE-QYQ-C-23×1600	2.3	1600
	ATE-QYQ-C-23×1800	2.3	1800
	ATE-QYQ-C-23×2300	2.3	2300
	Coated,with spike		
	ATE-QYQ-D-23×1600	2.3	1600
	ATE-QYQ-D-23×1800	2.3	1800
	ATE-QYQ-D-23×2300	2.3	2300
	Uncoated,alligator,without spike		
	ATE-QYQ-E-23×1600	2.3	1600
	ATE-QYQ-E-23×1800	2.3	1800
	ATE-QYQ-E-23×2300	2.3	2300
	Uncoated, alligator, with spike		
	ATE-QYQ-F-23x1600	2.3	1600
	ATE-QYQ-F-23x1800	2.3	1800
	ATE-QYQ-F-23x2300	2.3	2300
	Coated, alligator, without spike		
	ATE-QYQ-G-23×1600	2.3	1600
	ATE-QYQ-G-23×1800	2.3	1800
	ATE-QYQ-G-23×2300	2.3	2300
	Coated, alligator,with spike		
	ATE-QYQ-H-23×1600	2.3	1600
	ATE-QYQ-H-23×1800	2.3	1800
	ATE-QYQ-H-23×2300	2.3	2300

JIANGSU ATE MEDICAL TECHNOLOGY CO.,LTD

ADD: No.8,Lanxiang Road,Wujin Economic Development Zone
213000 Jiangsu China

E-mail: atemedical@xinhaigroup.com Web: www.ate-medical.com

 +86-519 8528 8988

 +86-519 8528 8060

A Xinwell Member Company

AteTec